REDAX®

DRENTECH™ EMOTRANS

IT

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DESCRIZIONE

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EN

DESCRIPTION

The Drentech® Emotrans systems are single-use thoracic drainage devices fitted with collection chambers, underwater seal valve, underwater columns for suction control and are predisposed for autologous blood recovery. The drainage systems can be used by suction or by gravity and are available in the adult and paediatric versions, with three chambers (single tube) and four chambers (double tube). The systems are composed of the following parts:

1. A collection chamber (1) divided into two graduated segments. The collection chamber has a maximum capacity of 2,150 cc and is divided into two sections that fill up in succession. The first section features a filter, with 200 µm porosity and large surface for the preliminary filtering of blood and macroaggregates elimination. In the double drainage version, the two sections of the collection chamber are independently graduated in order to allow separate loss evaluation.
2. A seal valve (2). This valve is mechanical, made by means of a silicone membrane and seal guaranteed when it is closed. Therefore, the "seal" is always guaranteed in any condition and there is no risk of blood contamination.
3. Two underwater columns for suction control (3) that allow regulating the negative pressure through the height of the water level inside (from -5 cmH₂O to -25 cmH₂O).
4. Back connector (8) for the connection to the Redax blood recovery bag, fitted with a connector with a pierceable membrane.
5. A label, situated on the side of the system, on which it is possible to put patient data and the starting time of autologous blood recovery.

INSTRUCTIONS FOR USE


The thoracic drainage systems are intended for:

- Evacuating air and/or fluids from the chest cavity or the mediastinum, preventing their accumulation;
- Re-establishing and maintaining normal pressure gradients in the thoracic cavity;
- Facilitating complete pulmonary expansion to restore normal respiratory dynamics;
- Collecting blood during postoperative phase through a preliminary filtering to eliminate microaggregates and surgical residuals.

The system is set for collected blood recovery, for immediate autotransfusion, thanks to the possibility of connection with the Redax blood recovery bag. Moreover, the system can be connected, by means of a suitable optional connection tube, to a separating device for the possible concentration and/or lavage of the erythrocytes before the reinfusion.

GENERAL WARNINGS

- The product is sterile if the package is intact. Do not use if the package has been opened or damaged.
- Dispose of the device after each single use; do not reuse.
- Avoid exposure to high temperatures and ultraviolet rays during storage.
- This device must be used by health care personnel trained for the use of thoracic drainage systems and blood recovery systems, with the supervision of a physician.

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WARNINGS/PRECAUTIONS

- The thoracic drainage system must be maintained in a vertical position, also during transport, and always below the fluid level of the pleural cavity (at least 50 cm). The system can be put in a horizontal position, as shown in figure 1, for a short time (ex. during patient transfer), without fluid blending between the different chambers or fluid release.
- The seal remains active also in horizontal position so it is not necessary to clamp the patient tube.
- In the event that the system overturns accidentally, the mechanical seal valve does not let the fluids container in the system blend with the collected blood. This latter remains unchanged and it is available for a subsequent reinfusion.
- Avoid kinks or squashing the patient drainage tube (9).
- Do not obstruct or cover the positive pressure release valve (5).
- Do not connect the system to the patient when the suction is already into operation.
- Thanks to its particular design, the system is stable and does not require a tube support for use on the floor. The system has a handle (10) to easier the transport of the drainage unit. The handle also allows hooking the system to the bed of the patient.
- Do not clamp the patient connection tube for a long time during the drainage and/or the transport.
- To use the system to collect and recover autologous blood, together with the Redax blood recovery bag, refer to the instructions for use of the latter.
- To use the Redax optional connection tube, to connect to a separating device, refer to the corresponding paragraph of this manual of instructions for use.
- As recovered blood is kept at room temperature, the time between the starting of the collection and the starting of the reinfusion must not be higher than 6 hours.

PREPARING THE DRAINAGE SYSTEM

The system is packed in double sterile wrapping. Remove the wrapping using an aseptic technique.

a. Preparing the seal valve

In the thoracic drainage systems, the seal valve is made by means of a silicone membrane so it does not require to be filled with water. This operation can be carried out, in an optional manner, to diagnose air loss (see the relevant paragraph in this instructions sheet). In this way, the valve guarantees air flow unidirectionality and the protection of the blood collected in the chamber in all situations.

b. Filling the suction regulating columns

1. Remove the protective cap (4) of the control columns (3).
2. Fill the control columns (3) with sterile water or saline solution up to the desired suction level (max. 25 cmH₂O). The water contained in the chamber will turn blue. The amount of liquid required depends on the desired pressure value; the approximate values corresponding to the pressures indicated are given below.

H₂O volume suction level:

10 cmH₂O 160 ml

15 cmH₂O 250 ml

20 cmH₂O 320 ml

25 cmH₂O 390 ml

3. Wait until the water level in both control columns (3) has aligned and corresponds to the desired suction level.
4. Close the cap (4) and close the front valve (11) as shown in the figure.
5. The water level in the control columns (3) determines the negative pressure values applied to the patient. The correct suction level must be chosen by a physician.

WARNING: it is important to check that the front valve (11) is correctly closed before connecting the suction source and starting to use the device. If the valve is not correctly closed, the columns would be unbalanced until partial fluid suction to the centralized source.

c. Connecting the drain tube

Connect the drain tube (9) to the thoracic catheter of the patient and then connect the other end to the collection chamber by means of the lock-in connector (12). In the version with four chambers, there are two connection tubes for connection to the thoracic catheters, identifiable by the numbering printed directly on the tube.

DEVICE OPERATION

a. Operation by suction

1. Connect a connection tube to the connector (6) to connect the drainage system to the suction source.

Note: The drainage system is equipped with an automatic "self-regulating" device (13) which allows obtaining moderate gurgling in the pressure regulating columns (3) in any vacuum condition available in hospital departments and without needing

any intervention by the staff. For this reason, the vacuum can directly be connected to the system without interposing other valves and/or regulators.

WARNINGS:

- Periodically check that suction in the control chamber (3) is continuous with moderate gurgling.
- Suction must be suspended in the following cases: while restoring the water levels or when you want to replace the system. After each suspension, check that suction is restored.
- Evaporation, especially at high temperatures, may cause the liquid level in the control columns (3) and in the water seal valve (2) to drop. This evaporation is normally negligible thanks to the "self-regulating" device. It is in any case advisable to periodically check the liquid level and top it up if necessary.

b. Operation by gravity

The drainage system can also be used by gravity by detaching the suction tube connected to the vacuum source from the connector (6). The automatic positive pressure relief valve (5) allows operation by gravity without requiring further operations.

c. Diagnostics of patient air loss

In the thoracic drainage systems, the seal is made by means of a silicone membrane and it does not require the filling of the valve with water. The diagnostics of possible air loss can be carried out through the filling of the lower chamber, as described below.

Note: This filling is optional and depends only on the physician need to evaluate these losses.

1. Fill a syringe with sterile water or physiological solution. It is not necessary to use the needle.
2. Insert the syringe in the connector (7) and fill the lower chamber, up to the broken line (about 45 ml). The water will turn blue thanks to a dye present in the chamber.
3. Remove the syringe; the connector (7) is provided with a valve that that stops access automatically when the syringe is removed.
4. Now it is possible to observe patient air loss.

WARNINGS

- The presence of air bubbles in the lower chamber of the seal valve (2) indicates an incomplete lung expansion. The reduction of the quantity of these bubbles in the course of time indicates a progressive lung expansion. When this expansion is complete, it will be possible to verify the absence of bubbles in the seal valve.
- Extensive bubbling in the seal valve (2) can be indicative of a serious pneumothorax or air loss in the connections. Carry out a careful checking and, if the phenomenon persists, call in a physician.

d. Reducing the suction level

The suction level can be reduced at any time if deemed necessary. The device allows easily and quickly reducing the suction level to half its original value. To carry out this operation, just open the protective cap (4). You will immediately notice that the characteristic bubbling disappears in the left-hand column while it remains unchanged in the right-hand column. That way, one of the two regulating columns is excluded thus reducing suction to about half the value originally set.

To restore the initial suction, just reposition the cap (4). The operations described can be carried out while the device is in normal operation and suction does not have to be suspended. Evaporation, especially at high temperatures, or considerable air leakage from the patient may cause the liquid level in the seal valve to drop. It is advisable to periodically check the liquid level and top it up if necessary.

e. Transporting the device

The device should preferably be transported in vertical position and below the level of the patient.

If necessary, the device can be placed in horizontal position as shown in Fig. 1. In this position, avoid mixing the liquids between the various chambers.

f. Collection and recovery of autologous blood

For bag connection operations and bag use, refer to the respective operating instructions.

WARNINGS

- As the recovered blood is kept at room temperature, the time intervening between the beginning of the collection and the beginning of the reinfusion must not exceed 6 hours: refer to the details shown on the label. In the event that this information is not reported, refer to the time of the end of the surgery.
- In the double-drainage version, only the blood in the first collection chamber can be recovered.
- Normally, it is not necessary to add coagulating substances. If the physician decides to use anticoagulant, on the basis of considerations concerning the clinical picture of the patient, it is

Quantità di Sangue - Blood quantity - Blutmenge - Quantité de sang - Cantidad de sangre - Hoeveelheid bloed - Quantidade de sangue - Blodmengde - Blodmängd - Verimäärä - Blodmængde - Ποσότητα αίματος - Blodmængde - Ilość krwi - Количество крови

Annotazioni - Remarks - Bemerkungen - Remarque - Observaciones - Opmerkingen - Anotações - Merknader - Anmärkingar - Huomautuksia - Bemærkninger - Σημειώσεις - Bemærkninger - Uwagi - Замечания

Identificazione paziente - Patient identification - Patientenidentifikation - Identification patient - Identificació del pacient - Identificatie patiënt - Identificação do paciente - Pasientidentifikasjon - Patientidentifying - Potilastiedot - Patient-id - Αναγνωρίση ασθενούς - Identifikation af patienten - Identifikacja pacjenta - Данные пациента

Ora inizio recupero - Recovery start time - Anfangszeit der Erholung - Temps de récupération - Tiempo inicio recuperación - Herstart start tijd - Inicio da recuperação - Klokkeslett for start av oppsamling - Tidpunkt för start av uppsamling - Talleenoton aloitusaika - Tidspunkt for start af prøvetagning - Ωρα αρχής ανάκτησης - Starttidspunkt for indsamlingen - Godzina rozpoczęcia przelaczania - Время начала восстановления

advisable to use ACD-A (in a dosage of 40-60 ml for each unit of collected blood or in an equivalent ratio of 1:7) because it is better metabolised. In patients with preserved hepatic functionality, the ACD-A anticoagulant solution does not require specific antagonists if the dosage of 150-180 ml is not exceeded. The calcium solution administration antagonises the effects of a possible excessive anticoagulant administration.

g. Replacing the unit

Replace the drainage system when the capacity of the collection chamber has been used up. To carry out this operation:

1. Prepare the new unit following the steps above so that it is ready for use.
2. Suspend suction and close the patient tube by means of the sliding clamp provided.
3. Disconnect the patient tube from the unit in use by means of the lock-in connector and replace it with the new one.
4. Re-open the clamp on the patient tube and restore suction.

DISPOSAL

Adopt adequate precautions for elimination and disposal of the device and comply with the provisions of the laws in force on biologically hazardous waste.

CONTRAINDICATIONS

No known contraindications exist for drainage use. For the recovery and reinfusion of autologous blood, see the Redax blood recovery bag operating instructions. In general, it is best to follow the "Universal precautions relating to blood and body fluids".

SINGLE-USE PRODUCT

Reuse may lead to alteration of performance and risks of cross-contamination.

STERILE - Ethylene oxide sterilized.

MATERIALS USED

Polystyrene (PS), Polypropylene (PP), Elastomer (S-E/B-S), Polyethylene (PEHD), ABS, PVC, LATEX-FREE MEDICAL DEVICE. PVC-FREE PACKAGE.

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