
DRENTECH™ mobile  
REDLINE

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### INSTRUCTIONS FOR USE

#### INTENDED USE OF THE DEVICE

The Drentech® Mobile REDLINE Thoracic Suction device (hereinafter referred to as "Vacuum Unit") is a portable unit with rechargeable batteries capable of generating adjustable suction. It is intended for use following thoracic surgery where a portable and independent vacuum unit is required.

#### DEVICE DESCRIPTION

The vacuum unit is a portable and adjustable suction system dedicated to use with compatible Redax thoracic drainage systems. The vacuum unit is composed of a plastic casing housing a user interface display (1), a keypad (2), a surface with fastening tabs (3) to secure the vacuum unit to the drainage system, a filter/valve connection port (4), a power supply connector (5) and a pin (6) to house the filter/valve pipette when it is not connected to the drainage system.

#### ACCESSORIES

The vacuum unit has the following accessories:

- A power supply (7) that allows recharging the batteries and continuous operation of the vacuum unit.
- A filter/valve (8) to allow safe connection to the thoracic drainage system (9).

The power supply is equipped with adapters for the various international electrical power supply standards; refer to the relative user manual for further details.

The vacuum unit may only be used in conjunction with a Redax single-use thoracic drainage system for which compatibility is expressly declared.

#### GENERAL WARNINGS/PRECAUTIONS

- **THE VACUUM UNIT IS SUPPLIED NOT STERILE AND MUST NOT BE SUBJECTED TO ANY STERILIZATION CYCLE.**
  - **The vacuum unit may only be used with Redax single-use thoracic drainage systems for which compatibility is expressly indicated. To this end, see the relative instructions for use.**
  - **Do not use the vacuum unit if its casing is clearly damaged.**
  - **Carefully follow the instructions to prepare the vacuum unit for use before using the system.**
  - **Do not expose the vacuum unit and the accessories to contact with liquids or atmospheric agents in conditions such as to exceed the IP44 degree of protection. For the cleaning operations, refer to the dedicated section below.**
  - **This device and each of its parts can be used in safe conditions only in the field of application and according to the procedures indicated in this instruction leaflet in correspondence to the type of product. The manufacturer declines all responsibility for improper use or different from that indicated.**
  - **The device may only be used by qualified and specialized medical staff (doctors and/or nurses), only in the field of application and according to the procedures indicated in this instruction leaflet in correspondence to the type of product. The manufacturer declines all responsibility for use by unauthorized persons.**
  - **The vacuum unit contains rechargeable batteries without memory effect, which can be recycled in all countries where differentiated waste collection is required and in any case in accordance with the local regulations in force.**
- NOTE:** to appropriately dispose of the batteries, they have to be removed from the container. The batteries may only be replaced/removed by REDAX personnel or in any case specialized persons authorized by REDAX.
- **During operation, do not cover the vacuum unit and/or the battery charger and/or the thoracic drainage system with clothes, gauze or other**

materials in order to prevent overheating.

- **The physician responsible and the assistant personnel must be aware of the possible complications related to high suction drainage for proper use of the vacuum unit with values higher than -30cmH<sub>2</sub>O (-3.0 Kpa) for which they are responsible in the last resort.**
- **If used at high altitude or on board aircraft, the performance of the vacuum unit might differ from the nominal.**
- **Do not expose the vacuum unit to high temperatures or fire because there is a risk of the batteries exploding.**
- **DO NOT DISCARD THE UNIT TOGETHER WITH THE BLOOD COLLECTION SYSTEM AFTER USE!**
- **To avert this eventuality, when the unit is turned off it flashes for a certain number of times at the same time emitting an acoustic alarm and showing the dedicated symbol (18) on the display.**
- **The unit has passed all the electromagnetic compatibility and electrical safety tests in compliance with the applicable European standards (CEI EN 60601-1-2): in any event, it is recommended to pay attention when the device is used in proximity of equipment that emits electromagnetic radiation during use.**

#### PREPARING FOR USE

##### *Preliminary battery charging*

The vacuum unit requires charging before it is used for the first time.

1. Take the vacuum unit out of the package.
2. Connect the power supply to the vacuum unit using the socket (5) on the rear of the unit.
3. Connect the power supply to a mains socket to start charging the batteries contained in the device.  
**Warnings for recharging and the batteries**  
The unit is equipped with rechargeable lithium polymer (LiPo) batteries. The vacuum unit is shipped with the batteries on low charge; it is advisable to put the device on charge until the batteries are fully charged. A complete charge from a completely flat battery takes about 1.5-2 hours.
4. Charging is indicated by the electric plug symbol (10) that appears on the display next to the battery charge status symbol (11) when the unit is on or by the battery symbol with an indication of the current charge status when the unit is off (12).

**NOTE:** The batteries will not be charged if they are damaged or the polarity is inverted or there is a short-circuit. In this event, a warning light will come on as described in the section "Warning lights and acoustic alarms".

##### *Preparing the vacuum unit*

Operate as follows to prepare the vacuum unit for use:

1. Prepare the Redax thoracic drainage system for use following the relative instructions for use and connect it to the drainage catheter.
2. Assemble the filter/valve unit (8) to the vacuum unit and connect the unit to the drainage system as shown in Figure 1. To this end, slide the drainage system on the upper surface of the vacuum unit taking care to fit the side and rear fastening tabs of the drainage system into place in the tabs on the vacuum unit.

#### USING THE DEVICE

##### *Using the vacuum unit*

Turn on the vacuum unit using the I/O button (15).

When the unit is turned on, it runs a functional test of the circuit with a solenoid valve excitation cycle (you should hear the valve opening and closing click) and pump activation for a few seconds. Should you not hear the solenoid valve opening and closing click and/or the noise of the running pump, the unit might be damaged. In this case, contact Redax Technical Service.

The battery charge (11), the suction value (17) and the suction unit of measure will be shown on the vacuum unit display. When the unit is turned on, it will be adjusted to "0 cmH<sub>2</sub>O", i.e. it will not perform any suction until the physician selects the desired negativity value. The unit is set to show the data in cmH<sub>2</sub>O; to show the data in kPa, contact Redax Technical Service.

Select the desired suction value by acting on the Arrow Up/Arrow Down buttons (14).

At this point, the vacuum unit will start operating intermittently, the frequency of activation depending on the amount of patient air and/or fluid leakage. You can change the suction value at any time by acting

on the respective selection buttons (14).

#### Recharging the batteries/power supply

When the batteries have discharged or whenever you want to restore full autonomy of the vacuum unit, you can recharge the batteries using the power supply (7). Connect the power supply to the mains socket and then plug it into the vacuum unit. Charging normally takes maximum 1.5-2 hours depending on whether or not you start with completely flat batteries. If the vacuum unit batteries are low and you want to continue drainage by suction, you can power the device by keeping it connected to the drainage system, therefore, during its normal use.

To disconnect the device from the power, disconnect the power supply.

#### Replacing the filter/valve unit

The vacuum unit protection filter/valve (8) must be replaced if necessary (e.g. prolonged use, fluids in the filter/valve). These accessories are supplied sterile and separate from the vacuum unit (see the catalogue or call a Redax representative). For replacement, disconnect the filter/valve in use and dispose of it and then connect a new filter/valve.

The above mentioned operations must be carried out with the vacuum unit off.

**NOTE:** The filtering membrane has a limited life and it is therefore advisable to replace it whenever you notice liquid and/or condensate on it, and in any case at least once a week.

The unit has been designed to prevent liquids from penetrating it. However, if you have any doubts about liquid having penetrated the unit, send it to Technical Service.

#### Operation by gravity

You can at any time suspend suction simply by turning off the vacuum unit using the dedicated button. At that point, drainage can continue by gravity without needing to disconnect the vacuum unit as it is equipped with an automatic positive pressure release valve which allows continuing drainage. It is however advisable to disconnect the pipette of the filter/valve in order to reduce pressure drops in air evacuation; the pipette can easily be housed on the pin (6).

At your discretion, suction can be restored at any time by turning the vacuum unit on again.

#### Operation at high suction

The vacuum unit normally operates in a range of suction values between -10 and -30 cmH<sub>2</sub>O (-1 to -3.0 KPa), universally recognized as "safety range" for thoracic applications. To increase suction, hold down the increase button for at least three seconds until the corresponding warning light comes on (see the relative section) accompanied by a short beep. At this point, the safety lock is released and you can increase the suction value up to a maximum of -60 cmH<sub>2</sub>O (-6 Kpa).

#### Keypad lock

After about 10 seconds of inactivity, the keypad will automatically be locked to prevent accidental and involuntary activation. To unlock the keypad, press the I/O button for 3 seconds until the padlock symbol (16) is shown on the display.

#### Turning off the device

To turn off the device, press and hold down the I/O button (15) for a few seconds until the display turns off.

**WARNING:** if the device is connected to the patient, before turning it off, make sure that it has expressly been requested that the patient be taken off drainage with a portable suction system.

#### WARNING LIGHTS AND ACOUSTIC ALARMS

The vacuum unit has 2 warning lights (see Fig. 1): the first (13) is located above the I/O button and the second is integrated in the display (i.e. backlit display). The unit also provides other information via the display.

Listed below are the various cases when one or more warning lights and/or symbols are shown.

**Power on:** the warning light (13) is GREEN and permanently on.

**Display:** in normal operating conditions, back-lighting is off. Pressing the I/O button (15), white back-lighting comes on, which makes it easier to read the display even in poor light conditions. In the event of an alarm (see below), the display lights up orange and flashes according to what is described below.

**Operation at high suction:** beyond -30 cmH<sub>2</sub>O, the warning light (13) is YELLOW and permanently on.

**Connection to power supply:** the electric plug symbol (10) is shown on the display.

**Low battery:** the display flashes orange twice every 30 seconds, the warning light

(13) is YELLOW and flashes and the battery symbol (11) flashes on the display. Each flash of the display is accompanied by a corresponding acoustic signal.

#### Circuit leak/suction fault/pump running continuously for more than 2 minutes:

the warning light (13) is YELLOW and flashes, the display flashes continuously, the pump continues running at maximum speed. This condition requires particular attention by the medical staff. It is advisable to promptly intervene if you notice this condition: if the pump inside the device remains active for a long period of time (2 hours), the temperature of the outer casing of the unit may exceed 41°C in some well localized areas (right-hand side).

POSSIBLE CAUSE	ACTION TO TAKE
<b>CIRCUIT LEAK</b>	Check all the system connections looking for possible partial or complete disconnections.
<b>VACUUM UNIT FAILURE</b>	Contact Redax Technical Service.
<b>HIGH PATIENT AIR LEAKAGE</b>	If the unit is on and connected to a patient with air leakage higher than 1,5-2 l/min, the pump may operate more or less continuously depending on the suction value selected.

**Error during recharging/power supply:** the warning light (13) is YELLOW and flashes, the display flashes YELLOW and the battery (11) and electric plug (10) symbols flash on the display. Each flash of the display is accompanied by a corresponding acoustic signal. When the unit is off, the battery and electric plug icons flash.

**Keypad lock:** when the keypad is locked, the padlock symbol (16) is shown on the display.

**DO NOT DISCARD warning:** when the unit is turned off, the warning light (13) flashes, the display flashes orange, the DO NOT DISCARD symbol is shown on the display and an acoustic signal is emitted at the same time.

**ACOUSTIC ALARM DISABLE:** the acoustic alarms can be disabled by simultaneously pressing the pressure regulating buttons until the crossed out speaker symbol appears on the display.

#### WHAT TO DO IF...

- **Liquid has penetrated the unit:** If any liquid has penetrated the unit, immediately suspend use and send it to Redax Technical Service.
- **The warning light (13) comes on and/or the display flashes YELLOW**

Check the nature of the fault referring to the information provided on the display.

#### In particular:

- If the battery symbol (11) flashes, it means that the batteries are running low and you need to recharge them as soon as possible otherwise the unit will automatically turn off. Connect the device to the mains as soon as possible using the power supply provided.
- If the electric plug symbol (10) flashes, it means that the power supply is connected but the battery is not properly being charged. Check the connection and functioning of the power supply and/or the condition of the batteries and replace them if necessary (in this case, contact Redax Technical Service).
- **The unit is accidentally dropped:** it might have been damaged, therefore, suspend use and send the unit to Redax Technical Service.

#### CLEANING AND MAINTENANCE

The unit does not require any routine maintenance by the user. Any operation must be carried out by authorized Redax technical staff. For cleaning and disinfection of the unit and its accessories, use products normally used in hospitals provided they are colourless, and follow the instructions below:

1. Moisten a cloth or gauze with a small amount of product.

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2. Wipe all the parts of the unit casing.

We suggest you use colourless skin disinfectants, solutions with a low chlorine content, methyl alcohol, pH neutral detergent solutions.

##### Warnings:

- Do not use Betadine or other disinfectant and/or detergent solutions that contain dyes to prevent damaging the casing or other parts of the unit and its accessories.
- Do not use acid or basic solutions other than those recommended to prevent damaging the exposed parts of the unit and its accessories.
- Always use detergent and disinfectant solutions with a cloth or gauze.

##### NEVER POUR THE PRODUCT DIRECTLY ONTO THE UNIT AND ITS ACCESSORIES.

##### DISPOSAL

As well as the constituent materials of the casing and the other support parts, the unit contains electronic components and rechargeable batteries. For this reason, the device must be disposed of in compliance with the regulations in force regarding disposal of electronic waste.

**The batteries may only be replaced/removed by REDAX personnel or in any case specialized persons authorized by REDAX.**


##### For the European Community

Dretech® Mobile Redline is marketed in compliance with the WEEE (Waste Electrical and Electronic Equipment) regulations.

The device is classified as APPLIED PART TYPE B in accordance with CEI EN 60601-1.

##### Symbols used

The unit has been labelled in compliance with the European regulations in force, in particular in accordance with CEI EN 60601-1 and correlated standards:

SYMBOL	DESCRIPTION
	Class II device (Ref. IEC 60417-5172)
	Applied part type B (Ref. IEC 60417-5840)
	Refer to the instruction manual (Ref. ISO 7010-M002)
	Dispose of waste electrical and electronic equipment separately (Ref. 2002/96/EC e 2006/66/CE)
	Direct current (Ref. IEC 60417-5031)
	Degree of casing protection (Ref. CEI EN 60529) - Protected against foreign bodies with a diameter greater than or equal to 1mm - Protected against water spray
	Date of manufacture (Ref. UNI CEI EN ISO 15223-1:2015)
	Serial number (Ref. UNI CEI EN ISO 15223-1:2015)
	Catalogue number (device code) (Ref. UNI CEI EN ISO 15223-1:2015)

In addition, the following symbols not specified by international directives or standards but universally recognized have been used:

SYMBOL	DESCRIPTION
	Universal recycling symbol
	Unofficial but universally recognized symbol referring to European Directive 2011/65/EU

##### TECHNICAL SPECIFICATIONS

###### VACUUM UNIT

Power: 12V 3S1P LiPo 3350 mAh  
Operating temperature: 15°C – 35°C  
Storage temperature: -20°C - +30°C  
Humidity (operating/storage): 30% - 70%  
Atmospheric pressure (operating/storage): 70.0 – 106.0 kPa  
Measuring tolerance: suction ±10%  
Weight: 800 g  
Degree of protection: Ip44

###### POWER SUPPLY

Model: Powerbox EXM30 5008  
Power: 15W  
Operating temperature: 0°C – 40°C  
Storage temperature: -40°C - +75°C  
Maximum power absorbed: 30W  
Maximum output current: 2.45A

###### LEGEND

- 01 Display
- 02 Keypad
- 03 Fastening seat
- 04 Filter/valve connection
- 05 Power supply connection
- 06 Pin for filter (when not used)
- 07 Power supply (accessory)
- 08 Filter/valve
- 09 REDAX thoracic drainage system
- 10 Plug symbol
- 11 Battery charge status symbol
- 12 Battery charging symbol when unit is off
- 13 LED
- 14 Up/down arrows
- 15 I/O button
- 16 Keypad lock padlock symbol
- 17 Set pressure value
- 18 DO NOT DISCARD symbol

Date of issue of latest version:  
See the last page: (REV.: XX-XXXX)

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