
Fig. 5

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#### DEVICE DESCRIPTION

The device called UNICO BASE is a percutaneous kit essentially made up of the following components:

1. A main body (1) with two accesses (2) with Luer-lock connection that allow connecting drainage devices such as syringes or collection bags according to the procedures described in this instruction leaflet. A selector (3) is fitted on the main body that can be turned 90° to allow drainage by gravity or injection of liquids such as rinsing solutions.
2. A one-way valve is integrated in the device body, which prevents fluid backflow to the patient in any condition.
3. A Verres needle (4) which is used for atraumatic introduction and fluid drainage. The above mentioned needle has a cannula with an atraumatic tip that protrudes from the needle end. This cannula can slide longitudinally until uncovering the cutting edge of the needle during introduction allowing easy penetration. Subsequently, the cannula returns to the extended position by means of a spring thus protecting the surrounding tissues from possible traumas. The needle is equipped with a grip (5) in which a visual indicator is fitted so that the user can see the position of the atraumatic protection stylet. A lock connector is fitted on the grip, which communicates with the cannula. The hollow cannula has a hole in the terminal section by means of which the fluids are directly aspirated both during positioning (check that the desired position has been reached) and drainage.

#### ACCESSORIES

The kits come standard with:

- A 2000 ml collection bag fitted with an air evacuation valve and an emptying valve (6).
- A syringe (7).

#### GENERAL WARNINGS AND PRECAUTIONS

- Carefully read the instructions for use before using the device.

- Use of the device is reserved for physicians and healthcare operators qualified to perform percutaneous access procedures and aware of the risks and possible implications of the above mentioned procedures.

- This device and each of its parts must be used in safe conditions only in the field of application and according to the procedures indicated in this instruction leaflet. The manufacturer declines all responsibility for any event deriving from improper use or different from that indicated.

- Check that the visual indicator of the protection stylet position slides smoothly and that it shows green when the stylet is in extended position. If the green indicator does not disappear during initial compression on the epidermis, stop introduction and replace the device.

- The device has been designed for use in association with Redax products and accessories. If different products and accessories are used, the physician is responsible for checking compatibility before use.

- The device may remain in place for maximum 29 days. The duration of the drainage procedure must be determined by the physician based on clinical considerations.

- Do not use the device if its wrapping or one of its parts is damaged.


- Patients diagnosed with coagulopathy or other blood diseases must be carefully examined for the risk of haemorrhages.

#### PREPARING FOR USE

Before using the system, carry out the preparation procedure as described below:

1. Check that the package is intact in order to ensure product functionality and sterility.
2. Open the packages contained in the kit using a sterile procedure.
3. Determine the needle insertion site, if necessary by means of X-ray or ultrasound scan.
4. Prepare a sterile field around the selected site taking the normal precautions for skin disinfection.
5. Take the Verres needle out of the package and check that the internal stylet slides smoothly.
6. Connect the main body of the device to the rear lock connector of

Fig.4b

EN

## UNICO BASE

Percutaneous kit

### INSTRUCTIONS FOR USE

#### GENERAL DESCRIPTION OF "CENTESIS"

The procedure consists of percutaneously inserting a needle/drainage catheter in order to draw out fluids (liquid and air) mainly but not exclusively from the pleural cavity of the patient. The purpose of this treatment is to relieve the symptomatology related to pneumothorax, including hypertensive, with pleural, pericardial and abdominal effusion and spillage of various nature; the procedure may also include drawing of fluids for diagnostic purposes or infusion of rinsing and/or therapeutic solutions.

#### INTENDED USE OF THE DEVICE

The UNICO system is prevalently intended for creation of an access, in a minimally invasive way, to the chest or abdomen of the patient in order to evacuate accumulated fluids and/or air. The device can be used for drainage by gravity or suction as well as for collection of fluid samples for diagnostic purposes.

The device can also be used for infusion of liquids such as saline rinsing solutions.

Fig. 3



Fig.4a

the needle and closed the protective caps (fig.1).  
7. Perform local anaesthesia of the site selected for needle introduction using the usual anaesthetics.

#### INTRODUCING THE NEEDLE

After doing the preparation described above, the needle can be introduced.

**WARNING:** Although the system is fitted with numerous safety devices, this operation is not without risks and it is therefore recommended to follow the procedure described below.

1. Firmly grip the device using both hands; use the hand on the needle grip to impart the insertion thrust. The other hand should preferably be positioned in proximity of the needle apex to allow maintaining the direction and avoiding excessive needle insertion after passing the dermis.
2. Rest the stylet tip on the patient's skin and exercise slight pressure until the green visual indicator disappears (fig. 2).
3. Very slowly continue with introduction until the green visual indicator reappears; this means that the end of the needle has penetrated and that the protection stylet has returned to its extended position performing its safety function (fig. 3) Needle entry and the consequent reappearance of the green indicator may be accompanied by a "click" due to the action of the spring that extends the protection stylet.

#### FLUID DRAINAGE: "PULL & PUSH" (suction)

Having positioned the needle, the fluid can be aspirated with the aid of a syringe and a collection bag. To connect the above mentioned accessories, open the corresponding caps on the two Luer-lock connectors on the main body.

Connect the syringe to one of the two connectors and the collection bag to the other one.

**NOTE:** The physician is free to choose which connector to use for the syringe or the bag, whatever is most comfortable.

Pull the syringe plunger (fig. 4a): the fluid will be aspirated into the syringe. Pushing on the plunger, the fluid will automatically be sent to the collection bag (fig. 4b) thanks to the one-way valves integrated in the device body and the collection bag. Repeating the pulling and pushing operation in sequence, the fluid will be evacuated and collected in the bag.

**NOTE:** No additional operation is necessary and the selector must remain in the initial position (the knob is oriented according to the needle axis as shown in fig. 3).

The bag can be replaced or the syringe or both accessories removed at any time based on the drainage result.

When removing the accessories, closed the protective caps to prevent possible contamination. As the device has an integrated one-way valve, there are no risks of backflow to the patient during the accessory connection and disconnection operations.

#### FLUID DRAINAGE: COLLECTION BAG (gravity)

After positioning the needle, you can proceed with fluid drainage by gravity. In this case, operate as follows:

1. Connect the collection bag to one of the Luer-lock connectors and make sure that the other connector is closed with its cap.
2. Appropriately position the collection system below the level of the patient.
3. Turn the selector on the main body up to the end of travel corresponding to a 90° angle with respect to the initial position (fig.5). At this point, the needle and the collection system are directly connected without interposition of any valve.

**WARNING:** Use only the collection bag included in the kit or Redax spare parts as they are equipped with a one-way valve that prevents fluid backflow to the patient.

4. If you want to replace the collection bag, for example, when it is full, turn the selector back to the initial position and then replace the bag. Once you have replaced the bag, again turn the selector by 90°.

#### CONTRAINDICATIONS

From examination of the medical literature, it results that percutaneous access to the chest may in very rare cases induce iatrogenic pneumothorax. This complication has been connected with handling errors in performing thoracentesis and using large-gauge needles and drains.

There are no specific contraindications for use of this device.

#### CLEANING AND MAINTENANCE

The device is single-use and hence does not require any cleaning or maintenance.

**SINGLE-USE DEVICE** - Reuse may lead to alterations in performance and cross-contamination risks.

#### DISPOSAL

After use, the device and all its accessories may be contaminated with

biological material, therefore, they must be disposed of as "biologically hazardous" waste in accordance with the regulations in force in the country where the device is disposed of.

Date of issue of the last version:  
see last page : (REV.: XX-XXXX)

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