REDAX®

IT

KIT PERCUTANEO UNICO

versioni: "STANDARD, XL e MULTI"

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include fluids extraction for diagnostic purposes or infusion of solutions (micronized talcum in sterile injectable solution, saline solution).

INTENDED USE OF THE DEVICE

The UNICO system is mainly intended for creating a minimally invasive access to the various cavities of the patient so that accumulated liquids and air can be removed. The device can be used for drainage by gravity or suction connected to collection systems as well as for collection of fluid samples for diagnostic purposes. The device can also be used for infusion of liquids such as saline solution or micronized talcum in sterile injectable solution.

DESCRIPTION OF THE DEVICE

The UNICO device is a percutaneous kit available in different versions and in different diameters.

The kit is essentially made up of the following components:

1. A main body with three accesses (1). The first access (2), located at the rear, is for inserting a Verres needle or an atraumatic stylet. The two other accesses (3) have Luer-lock connectors and can be used for connecting drainage devices such as syringes, drainage bags, suction systems, etc. following the procedures described on this instruction leaflet. A selector (4) is fitted on the main body that can be turned 90° to allow drainage by gravity or injection of liquids such as rinsing solutions. In the "standard" version, a drainage catheter (5) with holes and a depth marker is securely fastened to the main body. In the version "standard", the catheter is connected to the main body by means of a butterfly connector (6) which allows interposing an extension line (7).
2. A Verres needle (8) for introducing the catheter. This needle has a stylet with an atraumatic tip that protrudes from the end of the needle. This stylet can slide longitudinally until uncovering the cutting edge of the needle during introduction to then return to its extended position by means of a spring. The needle is equipped with a grip (9) in which a visual indicator is fitted so that the user can see the position of the atraumatic protection stylet.
3. The devices are available in the "standard" version where the catheter is straight and in the "multi-use" version where the catheter has a pig-tail end.

ACCESSORIES

The kits normally contain a 2000ml collection bag (11) fitted with an air evacuation valve (14) and an emptying valve (15).

On request, other accessories can be provided:

- Asyringe (10).
- A single-use scalpel (13).
- A drainage system connector (12).
- A single-use clamp (16), etc.

The extension line (7) is included in the respective product versions.

GENERAL WARNINGS AND PRECAUTIONS

- Carefully read the instructions for use before using the device.
- Use of this device is reserved for physicians and healthcare operators qualified to perform percutaneous access procedures and aware of the risks and possible implications of these procedures.
- This device and each of its parts must be used in safe conditions only in the field of application and according to the procedures indicated in this instruction leaflet. The manufacturer declines all responsibility for improper use or different from that indicated.
- Check that the visual indicator of the protection stylet position slides smoothly and that it shows green when the stylet is in extended position. If the green indicator does not disappear during initial compression on the external surface of the thoracic wall, do not continue introducing the needle and replace the device.
- The device has been designed for use with Redax products and accessories. If different products and accessories are used, the physician is responsible for checking compatibility before use.
- The device may remain in place for maximum 29 days. The duration of the drainage procedure must be determined by the physician based on clinical considerations.
- Do not use the device if its wrapping or one of its parts is damaged.
- Patients diagnosed with coagulopathy or other blood diseases must be carefully examined for the risk of haemorrhages.
- Do not in any way handle the drain or make extra holes with the aid of cutting tools.
- If securing the drain with a stitch, be extremely careful not to suture the tube or cut it with needles or sharp tools. Do not over-tighten the suture as this might cause the drain to break.
- After positioning the drain, check that the perforated section is fully inserted in the wound and that the connections are

EN

UNICO PERCUTANEOUS KIT

versions: "STANDARD, XL and MULTI"

INSTRUCTIONS FOR USE

GENERAL DESCRIPTION OF "CENTESIS"

The procedure involves percutaneous insertion of a drainage catheter in order to remove fluids (liquids and air). All the versions are dedicated mainly but not exclusively to use in the pleural cavity, abdomen, pericardium, etc. of the patient. The purpose of this treatment is to relieve the symptomatology related to - purely by way of example - pneumothorax including hypertensive, pericardial effusion and spillage or of different nature such as pleural or gastric; the procedure may also

- perfectly sealed.
- **Manually remove the drain pulling it gently avoiding abrupt manoeuvres. Do not use metal instruments, such as forceps or devices that may cause the drain to break.**
- **If the drain is left in place for very long periods of time, it may be difficult to remove. Be extremely careful during the removal operations.**
- **Once the drain has been removed, check that it is fully intact, an essential condition to exclude that drain fragments have been left in the cavity caused by accidental rupture.**
- **Single-use device. Reuse may lead to alterations in performance and cross-contamination risks.**

PREPARING FOR USE

Before using the system, carry out the preparation procedure as described below:

1. Check that the package is intact in order to ensure product functionality and sterility.
2. Open the package using a sterile procedure.
3. Determine the needle insertion site.
4. Prepare a sterile field around the selected site taking the normal precautions for skin disinfection.
5. Take the main body with attached catheter out of the package; the catheter is supplied assembled to the body also in the version with butterfly connectors.
6. In the version with pig-tail, slide the rigid tube (Fig. Pig-tail) present on the catheter towards the end of the catheter until it is completely straightened out.
7. Insert the needle through the rear access (2) until the grip is securely in place. The needle and its protection stylet must protrude from the catheter end by a few millimetres.
8. In the version with pig-tail, extract the rigid tube from the catheter.
9. Perform local anaesthesia of the site selected for needle introduction using the usual local anaesthetics.
10. If deemed necessary by the physician in order to facilitate introduction of the catheter, make a small incision in the skin.

WARNINGS

In cases of thoracentesis, the access point for needle insertion should be near the upper edge of the lower rib in order to avoid the neurovascular bundles located at the lower edge of each rib. For all other accesses, carefully check the position by means of radiographic or ultrasound examination before insertion.

NEEDLE INTRODUCTION AND CATHETER POSITIONING

After doing the preparation described above, the needle can be introduced.

WARNING: Although the system is fitted with numerous safety devices, this operation is not without risks and it is therefore recommended to follow the procedure described below.

1. Firmly grip the device with both hands. Use the hand on the main body to push in and guide the needle and hold its grip in contact with the device body. The other hand should be placed near the tip of the needle; gripping the catheter with the thumb and index finger, you can avoid inserting the needle too far into the patient's pleural cavity once you have gone through the skin.
2. Rest the needle tip on the patient's skin and exercise slight pressure until the green visual indicator disappears (Fig. 1).
3. Very slowly continue with introduction until the green visual indicator reappears; this means that the tip of the needle has entered the pleural space and that the protection stylet has returned to its extended position performing its safety function (Fig. 2). Entry into the cavity and consequent reappearance of the green indicator may be accompanied by a "click" due to the action of the spring that extends the protection stylet.
4. Move the device grip in order to orient the insertion direction of the catheter. To check that insertion has occurred in the desired site, i.e. in correspondence to the area where the fluid to be drained has accumulated, an exploratory sample can be taken. To this end, connect the syringe to the Luer-lock connector on the needle removing the corresponding cap, pull back the plunger and visually check (Fig. 3) that the fluid flows into the syringe; if so, remove the needle and position the catheter as shown in Figure 4.
5. Start removing the needle and at the same time completely position the catheter as shown in Fig. 4. In the version with pig-tail, when you remove the needle, the catheter will resume its rounded profile thanks to the shape memory of the material.

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Having positioned the catheter, the fluid can be drawn out of the pleural cavity with the aid of a syringe and a collection bag. To connect these accessories, open the corresponding caps located on the two Luer-lock connectors on the main body. Connect the syringe to one of the two connectors and the drainage bag to the other.

NOTE: The physician is free to choose which connector to use for the syringe or the bag, whatever is most comfortable. Pull the

syringe plunger (Fig. 5a): the fluid will be drawn out of the pleural cavity into the syringe. Pushing on the plunger, the fluid will automatically be sent to the collection bag (Fig. 5b). Repeating the pulling and pushing operation in sequence, the fluid will be evacuated and collected in the bag.

NOTE: No additional operation is necessary and the selector must remain in the initial position (the knob is oriented according to the catheter axis).

The bag can be replaced or the syringe or both accessories removed at any time based on the drainage result. When removing the accessories, it is recommended to close the protective caps in order to prevent possible contamination. In any event, should the caps remain open, there is no risk of backflow to the patient.

FITTING THE EXTENSION LINE (for the relative version)

The device version with butterfly connector allows fitting an extension line between the catheter and the main body. That way, the drain can easily be secured and kept in place for a long time providing greater comfort for the patient. To fit the extension line, clamp the catheter using the single-use clamp (16) included in the package. Then disconnect the butterfly connector and connect the line as shown in Fig. 6. Once you have completed the operation, the clamp can be opened and you can secure the catheter and connect the fluid collection accessories.

LIQUID DRAINAGE: COLLECTION BAG

After positioning and securing the catheter, you can proceed with liquid drainage by gravity. In this case, you need to connect the drainage bag included in the package to any one of the connectors available; the drained liquids are collected in the bag by effect of gravity.

To perform fluid drainage, operate as follows:

1. Connect the collection bag to one of the Luer-lock connectors and make sure that the other connector is closed with its cap.
2. Appropriately position the collection system below the level of the patient.
3. Turn the selector on the main body up to the end of travel corresponding to a 90° angle with respect to the initial position (Fig. 7). At this point, the catheter and the collection system are directly connected without interposition of any valve.

WARNING: Use only the collection bag included in the kit or Redax spare parts as they are equipped with a one-way valve that prevents fluid backflow to the patient.

4. If you want to replace the collection bag, for example, when it is full, turn the selector back to its initial position and then replace the bag. Once you have replaced the bag, again turn the selector by 90°.

FLUID DRAINAGE: DRAINAGE SYSTEM

After positioning and securing the catheter, you can proceed with fluid collection by gravity or suction using a dedicated postoperative collection system (e.g. silicone tank, thoracic drainage system, etc.).

In this case, you need to connect the collection system to any one of the connectors available using the dedicated connector included in the package (12).

To connect a postoperative drainage device, operate as follows:

1. If you want to use a thoracic drainage system, connect the connector included in the package to one of the Luer-lock connectors.
2. Prepare the drainage system following the instructions for use of the respective products and then connect the system to the free end of the connector.
3. Turn the selector on the main body up to the end of travel corresponding to a 90° angle with respect to the initial position. At this point, the catheter and the collection system are directly connected without interposition of any valve.

WARNING: Preferably use Redax products as they are equipped with a one-way valve that prevents fluid backflow to the patient. If connecting drainage devices of different manufacturers, make sure that they are equipped with a water or other type of one-way valve in order to ensure proper evacuation from the pleural cavity.

4. If you want to replace the collection system, for example, when it is full, turn the selector back to its initial position and then replace the bag or drain. Once you have completed the operation, again turn the selector by 90°.

INFUSION OF THERAPEUTIC SOLUTIONS

After positioning and securing the catheter, you can use it as access for infusion of solutions (micronized talcum in sterile injectable solution, saline solution).

To do this, operate as follows:

1. Connect the syringe containing the liquid to be infused to one of the Luer-lock connectors and make sure that the other connector is closed with its cap.
2. Turn the selector on the main body up to the end of travel

Fig. PNX

Fig. A


Fig. B



Fig. 7

corresponding to a 90° angle with respect to the initial position. At this point, the catheter and the syringe are directly connected without interposition of any valve and you can proceed with infusing the solution to the patient.

- Once you have completed the operation, turn the selector back to its initial position and disconnect the syringe. If you need to continue with drainage, you can carry out the operations described above (pull-push, drainage by gravity or suction); see above for details.

TREATMENT OF HYPERTENSIVE PNEUMOTHORAX

After doing the preparation as described above, the needle can be introduced.

WARNING: Although the system is fitted with numerous safety devices, this operation is not without risks and it is therefore recommended to follow the procedure described below.

- Firmly grip the device with both hands. Use the hand on the main body to push in and guide the needle and hold its grip in contact with the device body. The other hand should be placed near the tip of the needle; gripping the catheter with the thumb and index finger, you can avoid inserting the needle too far into the patient's pleural cavity once you have gone through the skin.
- Rest the needle tip on the patient's skin and exercise slight pressure until the green visual indicator disappears (PNX-Fig.A).
- Very slowly continue with introduction until the green visual indicator reappears; this means that the tip of the needle has entered the pleural space and that the protection stylet has returned to its extended position performing its safety function. Entry into the cavity and consequent reappearance of the green indicator may be accompanied by a "click" due to the action of the spring that extends the protection stylet.
- Move the device grip in order to orient the insertion direction of the catheter. To check that insertion has occurred in the desired site, i.e. in correspondence to the area where air has accumulated, and to perform a first chest decompression, open the Luer-lock on the rear of the needle (PNX-Fig.B).
- After checking correct positioning, start removing the needle and at the same time completely position the catheter (PNX-Fig.C). In the version with pig-tail, when you remove the needle, the catheter will resume its rounded profile thanks to the shape memory of the material.
- Connect the bag to one of the two Luer-lock connectors on the catheter body and turn the selector on the main body up to the end of travel corresponding to a 90° angle with respect to the initial position (PNX-Fig.D). At this point, the catheter and the collection system are directly connected without interposition of any valve.

WARNINGS:

- The cap (14) on the bag provided must remain OPEN to allow correct evacuation of air from the patient.
- Use only the Redax collection bag provided as it is equipped with a one-way valve that prevents fluid backflow to the patient as well as an air evacuation valve.

CONTRAINDICATIONS

From examination of medical literature, it results that the thoracentesis procedure may in very rare cases induce iatrogenic pneumothorax. This complication was related to manoeuvring errors when performing thoracentesis and the use of high gauge needles and drainage systems.

There are no specific contraindications for use of this device.

CLEANING AND MAINTENANCE

The device is single-use and hence does not require any cleaning or maintenance.

DISPOSAL

After use, the device and all its accessories may be contaminated with biological material, therefore, they must be disposed of as "biologically hazardous" waste in accordance with the regulations in force in the country where the device is disposed of.

Date of issue of latest version:
See the last page: (REV.: XX-XXXX)

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