


Fig. 5

### UNMOUNTING

CE 0123

Redax S.p.A.  
Via Galileo Galilei, 18  
Poggio Rusco (MN) Italy

## IT

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## **PORTABLE THORACIC VACUUM UNIT WITH AIR LEAKAGE AND INTRAPLEURAL PRESSURE DISPLAY**

### **INSTRUCTIONS FOR USE**

#### **GENERAL DESCRIPTION**

The Dretech® Palm EVO device (hereinafter referred to as "Unit") is a portable vacuum unit with rechargeable batteries, and is equipped with a display to show the following information:

- Indication of the "instantaneous" air leakage corresponding to the last minute of operation.
- Mean air leakage value of the patient in the last hour of measurement.
- Indication of the hours of device operation, calculated starting from the moment of activation.
- History of the air leakages and minimum and maximum intrapleural pressures relating to the last 99 hours of operation and available in numerical and graphic form.
- Measurement and storage of the minimum and maximum intrapleural pressure values of the patient relating to the last minute of operation.
- Real-time measurement of the minimum and maximum intrapleural pressure values for each respiratory act of the patient.

The Unit is also capable of generating suction that can be adjusted independent of central vacuum and/or energy sources. This function is used in all cases where a vacuum source is required and you are not connected to the central suction system or it is difficult to access and/or adjust.

The Unit can be powered and/or the batteries recharged using the power supply provided.

The Unit can transfer data to a USB storage device in order to view and/or store them on a PC. The USB storage device can be connected to the dedicated port on the Unit using the adapter provided, and the data can be viewed with standard spreadsheet programs or the RedaxPlot software provided.

#### **INTENDED USE**

The Unit is mainly intended for use in thoracic drainage procedures following traumas and/or surgery to create controlled suction and for objective measurement of the air leakages and intrapleural pressure of the patient. The Unit must be used in conjunction with compatible Redax collection systems and with its own sterilizing filter. Being independent of vacuum and energy sources, the Unit is particularly suitable in all cases where full patient mobility needs to be maintained.

**NOTE:** The device provides a measurement of the main parameters to help physicians evaluate the

EN

postoperative course of the patient. These parameters are represented in numerical or graphic form with a certain degree of approximation (see "Technical Specifications") and facilitate but do not under any circumstances substitute the evaluation of the attending physician and his or her decisions in managing the postoperative course.

#### DEVICE DESCRIPTION

The Unit (Figure 1) is made up of the following parts:

1. A keypad to select the basic functions of the device (01).
2. A display to view the basic functions (02).
3. A colour touch screen to view and select advanced functions (03).
4. Amini USB connector to transfer data to a PC (04).
5. A connector to connect the power supply (05).
6. A connector on the side to host the sterilizing filter (06).

#### ACCESSORIES

The Unit can be supplied individually or as part of a complete kit that may comprise the following accessories.

1. A power supply with multiple plugs and an extension line.
2. Batteries.
3. A mini USB-USB adapter
4. USB storage device (flash drive) containing the ancillary RedaxPlot software for data transfer from the memory of the Unit.

The Unit must be connected to a dedicated compatible Redax collection system in order to perform its function. The connection must be made using a dedicated protection filter (07) fitted with an automatic positive pressure release valve.

#### GENERAL WARNINGS AND PRECAUTIONS

- **The Unit may only be used in conjunction with Redax single-use collection systems for which compatibility is expressly indicated (see the relative instructions for use). If the device is not used in combination with compatible collection chambers, the performance could result partially or wholly compromised.**
- **Do not use the Unit if its casing or one of its parts is clearly damaged.**
- **Carefully read the instructions for use before using the Unit.**
- **Use the Unit in closed places, including any transport means (ambulance, helicopter ambulance, etc.). The device is not suitable for use in an open environment. Do not expose the Unit and its accessories to atmospheric agents.**
- **This device and all the accessories and/or optional devices must be used in safe conditions only in the field of application and according to the procedures indicated in this instruction leaflet. The manufacturer declines all responsibility for any event deriving from improper use or different from that indicated.**
- **The Unit may only be used by qualified medical and nursing staff that are aware of the possible implications related to thoracic drainage operations.**
- **This device is not suitable for use at high altitudes or in hyperbaric chambers; for these applications additional calibration might be necessary.**
- **Do not expose the device to high temperatures or fire because of the risk of the batteries exploding.**

#### USING THE DEVICE FOR THE FIRST TIME

When using the Unit for the first time, perform the preliminary operations described below.

1. Unpack the Unit.
2. Fit the batteries contained in the kit in their housing as shown in the illustrations 1-2-3.
3. Check the type of mains outlet and insert the corresponding plug.
4. Connect the power supply to the mains.
5. Connect the power supply to the vacuum unit to start charging

the batteries contained in the device. Charging is indicated by the electric plug symbol (12) that appears on the display and the battery symbol (11) whose segments appear sequentially.

#### WARNINGS

- The batteries will not be charged if they are damaged or the polarity is inverted or there is a short-circuit. In this event, a warning light will come on as described in the section "Warning lights and acoustic alarms".
- Before using the device for the first time, charge the batteries for at least 12 hours. If the Unit has been stored for a long period of time or has not been used for more than 6 months, the batteries must be charged for at least 24 hours in order to restore their useful capacity. In any event, when recharging batteries that have not been used for a long time, the display might indicate that they are fully charged long before they actually are. This is perfectly normal and depends on the initial condition of the batteries. Therefore, ignore the indication and continue charging the batteries for the times indicated above. The batteries will be restored to their full capacity after 2-3 charge and discharge cycles during normal use of the device.
- The power supply is designed to adjust to all mains voltages and frequencies, nevertheless it is advisable to check compatibility before connecting and to use an appropriate plug. In case of doubt, do not connect the power supply.
- Under no circumstances, especially when charging the batteries, let the power supply come into contact with flammable liquids or substances, and do not handle it with wet hands or cover it with cloths, gauze, clothing or other.
- In the event of faults when plugging in to the mains or while the batteries are being charged, immediately unplug and consult the section "Troubleshooting" in this instruction leaflet.

#### PREPARING FOR USE

- Connect the protection filter (07), first fitting it in the connector on the Unit and then fitting the other end in the dedicated connector on the collection system.
- Place the Unit on the preselected Redax single-use collection system. Use the hole underneath the Unit to fasten it on the dedicated pin on the collection system.
- Firmly press until you hear a click. To remove the unit, follow the instructions in the relative paragraph.
- When done, the system fastened correctly and ready for use should look as shown in Fig. 2.

#### TURNING THE UNIT ON AND OFF

The system is turned on and off by pressing the ON/OFF button (08).

At power on, the Unit runs a diagnostic test during which:

1. All the warning lights come on simultaneously.
2. A short beep is sounded.
3. The display and touch screen come on (they should look as shown in Fig. 3A).
4. Once the device completes the switch on procedure, the light indicator, placed above the switch on button, will be green.

#### RESETTING THE UNIT [New patient]

The Unit can store the mean air leakage and intrapleural pressure values of the patient in numerical and graphic form for a total of 99 hours of operation. When turning on the device, if you want to start counting and storage from the beginning, you need to clear the memory. To do this, operate as follows:

1. Unlock the device by pressing the ON/OFF button (08) holding it down for at least 3 seconds.
2. Access the Settings menu by tapping on the spanner symbol (20) on the touch screen.
3. Tap on the bin symbol (25); you will be asked to confirm the selection. At this point, you can still cancel the operation by tapping on the red cancel symbol (27) at the bottom left of the screen.
4. Confirming the operation by tapping the green confirmation symbol (28) at the bottom right, the Unit will completely be reset and ready for use on a new patient.

#### WARNINGS

- Each time you turn on the Unit, always check if any data from previous sessions is stored. The attending physician should evaluate whether this data is still valid, and if yes, new measurements can be made starting from the last hour. If not, the Unit must be reset.

#### USING THE DEVICE - BASIC FUNCTIONS

##### Unlocking the device

After about 10 seconds of inactivity, the keypad and the Settings

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menu on the touch screen will automatically be locked to prevent accidental and involuntary activation. To unlock the device, press the ON/OFF button (08) for 3 seconds until the padlock symbol (14) is shown on the display.

##### Selecting the suction value

- When the Unit is turned on, it always sets up for operation by gravity and "0" will therefore be shown on the display.
- Select the suction value by acting on the respective selection buttons (10A and 10B).
- At this point, the vacuum unit will start operating intermittently, the frequency of activation depending on the amount of patient air and/or fluid leakage. You can change the suction value at any time by acting on the respective selection buttons (10A and 10B).

##### Recharging the batteries/power supply

When the batteries have discharged or whenever you want to restore full autonomy of the vacuum unit, you can recharge the batteries using the power supply operating as described in the section "Using the device for the first time".

##### Powering the Unit during use

The Unit is powered by rechargeable batteries so that it can be used fully independently and hence allows the patient to move and walk about. The batteries periodically need to be recharged (see the technical data of the Unit) as described above.

If necessary you can connect the power supply to the Unit while it is being used. In this case, the power supply will recharge the batteries and at the same time power the Unit. In this way, the autonomy of the Unit is extended indefinitely.

**Warnings:** When the Unit is in use, part of the energy supplied is used to power the Unit, therefore, when the Unit is off, the battery life may be longer than declared.

##### Replacing the filter/valve unit

The vacuum unit protection filter (07) must be replaced if necessary (e.g. prolonged use, fluids in the filter, etc.). The filter/valve unit is supplied sterile and separate from the vacuum unit (see the catalogue or call a Redax representative). For replacement, disconnect the filter in use and dispose of it and then connect a new filter. The above mentioned operations must be carried out with the vacuum unit off.

**NOTE:** The filtering membrane has a limited life and it is therefore advisable to replace it whenever you notice liquid and/or condensate on it, and in any case at least once a week.

##### Operation by gravity

The Unit can operate by both suction and gravity. The default setting is operation by gravity and when the device is turned on "0" will be shown on the display. In standard operating mode, the Unit can measure both air leakages and the intrapleural pressure of the patient. The intrapleural pressure may slightly be affected by suction pump activations, necessary to measure air leakages, especially on the maximum values.

If you do not want the suction pump to activate, operate as described in the section "Advanced functions".

##### Warnings:

- It is always up to the attending physician to decide whether to perform drainage by suction or by gravity.
- If you select different operating modes or different suction levels, always make sure that the corresponding warning light is on.
- If the Unit is on but not connected to any accessory, operation by gravity may give rise to intermittent activation of the suction pump. This is quite normal, as the system is designed to operate connected to the filter and to a collection system and it should therefore not be considered a malfunction.

##### Operation at high suction

The vacuum unit normally operates in a range of suction values between 0 and -30cmH<sub>2</sub>O (0 to -3.0 KPa), universally recognised as "safety range" for thoracic applications. When this maximum value is reached, a safety lock is activated preventing any further involuntary or accidental increase. If the attending physician decides to increase the suction value to over -30cmH<sub>2</sub>O, hold the pressure increase button (10A) pressed down for at least 3 seconds to deactivate the safety lock. You will hear a short beep and the power light (29) will turn yellow. At this point, you can increase the suction value to a maximum of -60cmH<sub>2</sub>O (-6 KPa).

##### Symbols, warning lights and acoustic alarms

There are two warning lights on the vacuum unit (Fig. 3A) situated on

the left of the display: one (29) at the bottom in proximity of the ON/OFF button and the other (30) above it. The Unit also provides other information via the LCD display and the touch screen.

Listed below are the various cases when one or more warning lights and/or symbols are shown.

**Power on:** The warning light (29) is GREEN and permanently on.

**High suction mode:** The warning light (29) is YELLOW and permanently on.

**Pump OFF mode:** The light indicator [29] is lightning yellow and fixed. ON the touch screen will appear, both on home page and in the other pages, the OFF sign.

**Battery charge/power supply OK:** The electric plug symbol (12) is shown on the display.

**Low battery:** The warning light (30) is YELLOW and flashes and the battery symbol (11) flashes on the display.

**Circuit leakage:** The warning light (30) is YELLOW and flashes and the word LEAKAGE (13) is shown on the display. This indicates very high air leakage or a possible disconnection in the circuit connected to the patient.

**Battery charge/power supply not OK:** The warning light (30) is YELLOW and flashes and the battery symbol (11) and the electric plug symbol (12) flash on the display.

**Pressure alarm:** If the intrapleural pressure remains constant at a certain value for more than 11 consecutive minutes, the yellow warning light (30) flashes as does the bar indicator (18) on the touchscreen.

This alarm may indicate:

- Occlusion of the patient catheter or tube
- Occlusion of the filter/valve
- Vacuum unit incorrectly connected to the drain.

**Keypad lock:** When the keypad is locked, the padlock symbol (14) is shown on the display.

For the meaning of the symbols and the indications shown on the LCD display and the touchscreen, respectively, refer to Fig.3A and Fig. 3B and to the legend. You can also set an acoustic alarm connected with activation of the warning light (30). For this operation, see the section "Settings menu".

## USING THE DEVICE - ADVANCED FUNCTIONS

### AIR LEAKAGE MEASUREMENT

This function allows objective viewing of the air leakages of the patient in real-time, the mean hourly value, and the stored historical data of the postoperative course. The real-time leakages are always shown in the centre of the display expressed in ml/min. This data is updated every minute and stored in the Unit. You can also view the historical trend of the postoperative course; this is given as mean hourly value and is updated hour by hour and stored in the Unit. The historical trend of the air leakages can also be viewed in graphic form directly on the screen.

#### Viewing the leakages in real-time

This information is displayed on the touch screen (03) and you can go from this view to the mean hourly value view simply by tapping the centre of the screen near the numerical value. The first time the Unit is turned on and after each reset operation, the display will show the following information:

- Air leakage indicator in the centre of the screen, which will show "0.0l/min".
- Operating hour counter at the top left, which will show "0". If values different from the above are shown, it means that some data have remained stored from previous use. If these data relate to the same patient and/or you want to keep them, the Unit will continue measurement and storage from these data. If you want to clear the stored history, operate as described in the section "Resetting the unit".

The Unit will start measuring the air leakages and every minute of operation display the leakage value expressed in ml/min. This data is updated minute by minute.

### Viewing the mean hourly leakage value

This information is displayed on the touch screen (03) and you can go from this view to the real-time view simply by tapping in the centre of the screen. The mean hourly value is expressed in ml/min and is indicated by the green caption "AVG 1H". This value is updated every hour of operation.

### Viewing the air leakage history - numerical values

Every hour of operation the Unit stores the mean hourly air leakage value of the patient. It can store up to 99 hours, being the maximum capacity of the memory. The stored history can be viewed starting from the end of the first hour. To view the air leakage history in numerical form, tap the symbol at the bottom left of the Home page (21). Left and right arrows (23) will appear, which you can use to scroll the list of values stored in either direction. The hour counter will always show the measurement time relating to the value on the screen at that moment. To exit and go back to the Home page, tap the relative symbol (22) at the bottom right of the screen.

The system repeats the measurement, display and storage procedure for a maximum of 99 hours (about 4 days) after which no further values will be stored until the Unit is reset. This will be indicated by the operating hour counter which will show the value 99 flashing.

### Warning:

- It is advisable to periodically check the operating hour counter to see if the maximum memory capacity has been reached so as not to lose significant data.

### Viewing the air leakage history - graphic mode

The historical trend of the air leakages can also be viewed in graphic form in relation to the elapsed hours from application of the Unit. To access the graph from the Home page, tap the LEAK symbol (19) at the bottom right of the touch screen.

As shown in Fig. 4, the graph shows the air leakage values (in ml/min) on the VIOLET ordinates (left axis), the suction value set in the time interval considered (cmH<sub>2</sub>O or kPa) on the GREEN ordinates (right axis), and the time (in hours) on the bottom axis. The air leakage scale automatically adjusts to the maximum value measured at any time during the postoperative course, so that the physician can easily see the leakage trend over time.

Each page displays 8 hours of history. If the Unit contains data relating to a longer time interval, you can scroll the graph in both directions by tapping the relative arrows (23) at the bottom of the screen (Fig. 4).

The air leakage graph is plotted in correspondence to the suction value set, so that the physician can easily associate the respective values in order to obtain information on the postoperative course.

To go back to the Home page, tap the relative symbol (22) at the bottom right of the screen.

### INTRAPLEURAL PRESSURE MEASUREMENT

This function allows objective viewing of the intrapleural pressure values measured in the collection chamber of the drainage device. The intrapleural pressure is displayed highlighting the maximum and minimum values, typically corresponding to the patient's inspiration and expiration phases; this data is moreover shown in real-time and as historical values.

### Viewing the intrapleural pressure in real-time

On the touch screen on the Home page (Fig. 3A) there is a graphic indicator (18) resembling a luminous scroll bar. The two-colour band oscillates synchronously with the respiratory trend of the patient and displays pressures in a range between +15cmH<sub>2</sub>O (+1.5kPa) and -30cmH<sub>2</sub>O (-3kPa). That way, the physician can in real-time monitor the intrapleural pressure trend for each single respiratory act.

### Viewing the intrapleural pressure history - graphic mode

The historical trend of the intrapleural pressure is only available in graphic form as it is stored every minute.

To access the graph from the Home page, tap the LEAK symbol (19) at the bottom right of the touch screen (Fig. 3A). Starting from the air leakage graph, tap the PRESS symbol (24) at the bottom left of the touch screen.

The graph shows the pressure values (in cmH<sub>2</sub>O or kPa) on the vertical axis and the time (in hours) on the horizontal axis (Fig. 5). The pressure scale is always in a viewing range between +15 and -30 cmH<sub>2</sub>O. In particular, the maximum pressure trend is represented by a RED line, while the minimum pressure trend by a YELLOW line.

Each page displays 8 hours of history. If the Unit contains data relating to a longer time interval, you can scroll the graph in both directions by tapping the relative arrows (23) at the bottom of the

screen.

To go back to the Home page, tap the relative symbol (22) at the bottom right of the screen.

### SETTINGS MENU

You can make a few settings by accessing the dedicated menu. To access the Settings menu from the Home page, unlock the keypad with the ON/OFF button and tap the spanner symbol at the bottom centre of the touch screen (20).

The Settings menu has the following options:

#### - Delete stored data

See the section "Resetting the Unit".

#### - Turn off the suction pump [PUMP OFF]

If you want to perform drainage and/or evaluate the intrapleural pressure completely excluding any suction pump activation, you can turn it off. From the Settings menu, tap the PUMP ON/OFF button (26) and the pump will be deactivated. The word OFF will be shown at the top of the touch screen (31) and the yellow light indicator [29] so that you can see that it is off. You can always reactivate the pump by repeating the steps above and tapping the same button in the Settings menu.

**IMPORTANT:** When the suction pump is deactivated, you cannot make air leakage measurements.

#### - Activate acoustic alarms

The Palm unit is delivered with the acoustic alarms activated. You can de-activate them at any time by tapping the speaker button (31) on the touch screen. Activation of an acoustic alarm also makes the warning light (30) come on indicating that operator intervention is required (see the section "Warning lights and acoustic alarms"). The details of what the operator needs to do can be deduced from the indications on the display.

### TRANSFERRING DATA VIA THE USB STORAGE DEVICE

You can transfer the stored data to a USB storage device (flash drive) via the mini USB port (04) on the front of the Palm unit. You can use the flash drive provided or a similar one to carry out this operation.

1. Remove the screw that secures the protective cap of the mini USB port and then open the port by sliding the cap sideways.
2. Insert the adapter provided. Alternatively, you can use a similar adapter as long as it is the same type as the one provided.
3. Insert the USB flash drive in the adapter.
4. From the Home page, select to display one of the graphs available by acting on the respective buttons.
5. You will see a button appear bearing the harmonised USB symbol in red. Tapping this button, the data will be transferred from the Palm unit memory to the USB flash drive.
6. When transfer is complete, the USB symbol changes to green.
7. At this point, you can remove the USB flash drive and the adapter and close the protective cap.

The file transferred to the USB flash drive is type .csv and can be opened with a common spreadsheet program (e.g. Excel, OpenOffice-Calc, etc.).

### NOTE:

To correctly open the file, access the spreadsheet program and use the "File Open" option. Generally, the structure is automatically recognised and the file correctly opened by the most common software programs. Should you have difficulties, select only the separator ";" from the options available and nothing else.

To directly view the graphs, you can use the dedicated RedaxPlot software provided on the USB storage device or request it by sending a mail to the following address: info@redax.it. To use the RedaxPlot software, refer to the instructions for use.

### WHAT TO DO IF...

#### - Liquid has penetrated the unit

If liquid penetrates the unit, it will be expelled through a dedicated outlet at the bottom of the unit. If any liquid has penetrated the unit, immediately suspend use and send it to Technical Service and/or consult the maintenance manual.

#### - The warning light (30) comes on.

Check the nature of the fault referring to the information provided on the LCD display.

In particular:

- o If the battery symbol flashes, it means that the batteries are running low and you need to recharge them as soon as possible otherwise the unit will automatically turn off.

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- o If the electric plug symbol flashes, it means that the power supply is connected but the battery is not properly being charged. Check the connection and functioning of the power supply and/or the condition of the batteries and replace them if necessary.
- o If the LEAKAGE symbol appears, it means that the unit was unable to reach the suction value set because of a high air flow. This may be a sign of a disconnection in the patient circuit and you therefore need to check the connections and the actual air leakage of the patient.
- o If the bar indicator (18) on the touchscreen flashes, it means that the pressure measured has remained constant for more than 12 minutes. There may be various reasons for this:
  - The patient has no air leakage
  - The patient tube and/or catheter are occluded
  - The unit is incorrectly positioned on the drainage system

Do the necessary checks to identify the problem and take the appropriate action.

#### **- The unit is accidentally dropped**

It might have been damaged, therefore, suspend use and send it to Technical Service and/or consult the maintenance manual.

#### **- Not information is transferred to the USB key**

Remove the key and repeat the operation described in the respective paragraph.

#### CLEANING AND MAINTENANCE

The Unit does not require maintenance by the user. Any operation must be carried out by authorised Redax technical staff. For cleaning and disinfection of the unit and its accessories, use products normally used in hospitals provided they are colourless, and follow the instructions below:

1. Moisten a cloth or gauze with a small amount of product.
2. Wipe all the parts of the Unit casing.

We suggest you use colourless skin disinfectants, solutions with a low chlorine content, methyl alcohol, pH neutral detergent solutions.

#### Warnings:

- Do not use Betadine or other disinfectant and/or detergent solutions that contain dyes to prevent damaging the casing or other parts of the Unit and its accessories.
- Do not use acid or basic solutions other than those recommended to prevent damaging the exposed parts of the Unit and its accessories.
- Always use detergent and disinfectant solutions with a cloth or gauze.

#### NEVER POUR THE PRODUCT DIRECTLY ONTO THE UNIT AND ITS ACCESSORIES.

#### REMOVING THE VACUUM UNIT

To remove EVO from the drainage system, operate as follows:

1. Turn off the device.
2. Disconnect the side filter.
3. Grip the unit from the rear and extract it from its seat making a rotary movement as shown in Fig. 6.

#### CONTRAINDICATIONS

No contraindications for use of this device are known at the time of printing this instruction leaflet. Ultimately, the decision on which drainage procedure to use is up to the attending physician.

#### ELECTROMAGNETIC INTERFERENCE WITH OTHER DEVICES

The Unit has been subjected to all the tests as required by the regulations in force regarding Electromagnetic Interference and Susceptibility. Therefore, during normal use, it is not subject to influences nor interferes with other equipment in typical hospital environments.

#### Warnings:

- Do not use the Unit in proximity of Magnetic Resonance equipment as it contains metal parts that may interfere with operation of such equipment.

#### DISPOSAL

As well as the constituent materials of the casing and the other support parts, the unit contains electronic components and rechargeable batteries. For this reason, the device must be disposed of in compliance with the regulations in force regarding disposal of electronic waste.

#### For the European Community

Drentech® Palm Evo is marketed in compliance with the WEEE

(Waste Electrical and Electronic Equipment) regulations.

#### TECHNICAL SPECIFICATIONS VACUUM UNIT

Power: 6V 5xAA NiMH 2500 mAh

Operating temperature: 10 – 35°C

Storage temperature: 0 – 60°C

Measuring tolerance: suction ± 10%, air leakage ± 10%

Weight: 540 g (excluding batteries)

#### POWER SUPPLY

Model: FRM015-S12-I

Power: 100-240V

Maximum current input: 0,7 A

Operating temperature: 10 – 35°C

Storage temperature: 0 – 60°C

Maximum power output: 15W

Maximum vacuum power: < 0,3W

Maximum output current: 1.25A

**NOTE:** the 12V out battery charger completely fits with the PalmEVO 6V input. Indeed, the higher voltage, provided by the battery charger is needed for a proper charging and functioning of the unit.

#### LEGEND:

- 01 Keypad
- 02 LCD display
- 03 Touch screen
- 04 Mini USB connector
- 05 Power supply connector
- 06 Protection filter connector
- 07 Protection filter
- 08 ON/OFF button
- 09 Touch screen on/off
- 10 Suction selection arrows
- 11 Battery symbol
- 12 Electric plug symbol
- 13 "Excessive air leakage" symbol
- 14 "Keypad lock" symbol
- 15 Operating hour counter
- 16 Air leakage indicator
- 17 Mean hourly value indication
- 18 Intrapleural pressure display scale
- 19 Go to air leakage graph
- 20 Settings
- 21 View air leakage history in numerical form
- 22 Go to Home page
- 23 Scroll arrows
- 24 Go to intrapleural pressure graph
- 25 Clear data
- 26 Suction pump on/off
- 27 Cancel clearing the memory
- 28 Confirm clearing the memory
- 29 Power light
- 30 General warning light
- 31 Speaker

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