
13

DRENTech™ mobile
REDLINE

CE 0123

Redax S.p.A.
Via Galileo Galilei, 18
Poggio Rusco (MN) Italy

REDAX

DRENTech™ SIMPLY
DRENTech™ mini SIMPLY

IT

ISTRUZIONI D'USO

DESCRIZIONE

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vedere ultima pagina : (REV.: XX-XXXX)

intrathoracic pressure gradients; facilitating complete pulmonary expansion to restore normal respiratory dynamics. The systems can perform their function by gravity or suction.

Note: This device can be used in safe conditions only in the field of application and according to the procedures indicated in this instruction leaflet in correspondence to the type of product. The manufacturer declines all responsibility for improper use or different from that indicated.

GENERAL WARNINGS

- The product is sterile if the package is intact. Do not use if the package has been opened or damaged.
- Single-use device: dispose of the device after each single use; do not reuse.
- Avoid exposure to high temperatures and ultraviolet rays during storage.

WARNINGS/PRECAUTIONS

- **The drainage system must be held in vertical position, including during transport, and always below the fluid level of the pleural cavity (at least 50cm).**
- **Avoid looping or squashing the patient drain tube.**
- **Do not obstruct or cover the water seal chamber filling connector (5).**
- **Operation by suction is enabled by connecting the system to the centralised suction source of the hospital provided that it is equipped with a thoracic regulator and under direct responsibility of the attending physician. Alternatively, it can be connected to the Drentech RedLine portable vacuum unit.**
- **The fluid contained in the collection chambers of the drainage system must be disposed of and may absolutely not be used for infusion.**

PREPARING AND USING THE DRENTech SIMPLY SYSTEM

The system is packed in double sterile wrapping. In order to use the product, take it out of its wrapping using a sterile technique.

a. Filling the water seal

1. Prepare a preferably 60ml syringe or a receptacle containing sterile water or saline solution.
2. Fill the water seal through the filling access (4) at least up to the solid line, which is the minimum reference level and corresponds to about 45ml (also see the indication on the front panel). The water will turn blue thanks to a dye contained in the chamber.
3. When the device drains by gravity, the water seal acts as pressure gauge. On the relative column (11), you can qualitatively observe the intrathoracic pressure oscillations synchronous with the patient's respiratory acts.
4. Air bubbles may appear in the chamber underneath the water seal: in normal conditions, these indicate air leakage from the pleural cavity, whose entity can be qualitatively measured (see "e. Patient air leakage diagnostics").

b. Connecting the patient tube

1. Close the patient tube (6) with a clamp or similar and connect it to the thoracic catheter of the patient by means of the connector (12). Then connect the patient tube to the collection receptacle by means of the lock-in connector (7). Fit the connector until you hear a click indicating that it is correctly positioned. At this point, the system is ready to operate by gravity.
2. To connect the device to the centralised suction source, first of all prepare a thoracic pressure regulator with an adjustment range within the safety parameters for the intended application (standard range used: -10 to -25 cmH₂O = -1 to -2.5 kPa; maximum suction: -60 cmH₂O = -6 kPa). After adjusting the suction to the desired value, detach the dedicated connector (5) from the handle by turning it, tightly insert it into the filling access (4) and then connect the suction. The device ready for use should look as shown in Fig. 1.


3. For operation by suction, prepare a RedLine portable vacuum unit following the instructions given in the relative instruction leaflet. Insert Drentech Simply into the dedicated seat on the unit and then connect the protection filter to the relative connectors and start using the device by controlled suction.

c. Disconnecting the patient tube

This operation can be carried out to replace the collection receptacle or to better arrange the patient tube, for example, based on the bed configuration.

1. Close the patient tube with a clamp or similar.
WARNING: Never disconnect the quick-coupling before having closed the tube with a clamp or similar.
2. Disconnect the lock-in connector (7) by pressing on the sides and at the same time pulling at the two ends. Then carry out the desired operations, for example, replace the collection receptacle.

Fig.1

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VÝROBEK NEOBSAHUJE LATEX.

OBAL VÝROBKU NEOBSAHUJE PVC

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- Finally, refit the connector into its seat and press until you hear a click.

d. Resetting high negative pressure

During operation by gravity, high negative pressures may be created in the collection receptacle (this does not occur when using the vacuum unit). These pressures spontaneously balance out over time. If you want to manually reset the pressure, temporarily disconnect the lock-in connector (7) of the patient tube following the procedure described in c. above.

e. Patient air leakage diagnostics

- Air bubbles in the lower chamber of the seal valve (2) are an indication of air leakage and, consequently, incomplete lung expansion. A reduction in the amount of bubbles over time indicates a gradual reduction of the leakage and better lung expansion. When lung expansion is complete, there will be no bubbles in the seal valve.
WARNING: If the phenomenon persists, ask a physician to intervene after carefully checking the connections.
NOTE: Evaporation, especially at high temperatures, or considerable air leakage from the patient may cause the liquid level in the seal valve to drop. It is advisable to periodically check the liquid level and top it up if necessary.

f. Transporting the device

The device must be transported in vertical position and below the level of the patient. In the event that the device is accidentally dropped, the Drentech Simply design prevents or strongly limits any fluid leakage allowing you to restore normal operation once the device is put back into vertical position.

WARNING: If the device is dropped, the water seal is compromised, therefore, put it back into vertical position as soon as possible and check the water level in the seal valve.

DISPOSAL

Adopt adequate precautions for elimination and disposal of the device and comply with the provisions of the laws in force on biologically hazardous waste.

Two caps (13) are fitted on the handle, which can be used to close the device accesses in order to facilitate disposal. The caps can easily be detached by turning them.

CONTRAINDICATIONS

There are no known contraindications.

SINGLE-USE DEVICE. Reuse may lead to alterations in performance and cross-contamination risks.

STERILE - Sterilized by ethylene oxide

MATERIALS USED

See the product technical data sheet.

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