REDAX®

IT

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INSTRUCTIONS FOR USE

DESCRIPTION

The SMART DRAIN drains with grooved round profile and coaxial internal lumen are made of pure biocompatible silicone and are composed of:

- A grooved tube section. This exclusive design allows effective drainage by effect of the capillarity principle and prevents tube occlusion following traction and/or torsion.
- This section has a round lumen, coaxial to the grooved section, which allows separate air evacuation from the thoracic cavity drawing it directly from the distal end of the tube through specifically created holes.
- A transparent silicone tube section, sealed and connected to the grooved section, which follows an atraumatic profile by means of which you can connect to closed gravity or suction drainage systems.
- The absolutely even diameter makes the device suitable for use in laparoscopy and thoracoscopy. The transition between the grooved and the round profile occurs rapidly in order to prevent the formation of small-diameter ducts that could be obstructed by clots. Depth markers, graduated from the end of the radiopaque end, facilitate drain positioning and removal.

CONTRAINDICATIONS

The drains are intended for short-term use in clinical conditions that require evacuation of air and/or fluids that have built up in natural and/or newly formed cavities following surgery. They are suitable for use in association with gravity or suction drainage systems.

Note:

These devices can be used in safe conditions only in the field of application and according to the procedures indicated in this instruction leaflet in correspondence to the type of product. The manufacturer declines all responsibility for improper use or different from that indicated.

SUGGESTIONS FOR USE

Positioning the drain.

- 1. Before positioning the drain, carefully choose the type, position and size required. These parameters depend on the type of operation, the physiology of the patient and the experience of the surgeon.
2. Using an aseptic technique, carefully take the sterile thoracic drain out of the protective wrapping.
3. Position the drain in the wound following an aseptic procedure based on the protocols and practices in use at the hospital.
4. For correct positioning, use the depth markers available on all the drains.
5. Secure the drain with a skin stitch or a plaster.
6. Connect the positioned tube to the preselected drainage system according to the procedures indicated in the relative instruction leaflet. If necessary you can use the connector included in the package (in the models where provided) to make the connection to the drainage system.

Removing the drain.

- 1. Remove the tape or suture stitch or any other type of fastening.
2. Disconnect the drainage system.
3. Remove the drain by slightly pulling it and immediately apply a dressing on the wound.

In thoracic surgery, the best drainage efficiency is obtained by connecting the drain tube to the Drentech Simple 3 system with double-lumen tube and air/liquid separation system (code 10186/10187). In fact, the air and liquid separation that occurs in the drain tube continues in the patient tube up to the collection receptacle. That way, air evacuation can occur without any obstacles and continuously and more effectively.

WARNINGS/PRECAUTIONS

- 1. Before inserting the drains, identify the insertion site and choose the correct size.
2. The choice of the drain size is based on the type of disease and surgery and the experience of the physician.
3. It is strictly prohibited to handle the drain and make additional holes with any instrument. Only the length of the holed terminal section may be reduced, provided that a clean cut is made without slashes or burrs and that the remaining holed section is not less than half the original length in order not to reduce the draining capacity.
4. Once extracted, check that the drain is fully intact, an essential condition to exclude that any drain fragments generated by accidental rupture have remained in the cavity.
5. It is recommended not to use the drain for a period of more than 29 days.
6. After positioning the device, check that the grooved section is fully included in the wound and that the connections are perfectly sealed.
7. Remove the device by slightly pulling it by hand avoiding sharp movements. Do not use metal instruments, such as forceps or other that may cause the device to break.
8. If the device is left in place for very long periods of time, it may be difficult to remove. Be extremely careful during the

removal operations.

- 9. Do not use roller clamps or other metal devices to milk the silicone tubes.
10. Do not use abrasive instruments, such as wads of cotton or gauze. Milking the silicone tubes using this type of instrument may cause damage to the drain surface reducing its mechanical resistance.
11. Do not use any alcoholic or aggressive liquid to lubricate the tube surface before milking. Redax strongly suggests using a cream or Vaseline oil or any other non-aggressive lubricant to help milking the tube.

CONTRAINDICATIONS

There are no known contraindications.

GENERAL WARNINGS

The product is sterile if the package is intact. Single-use device. Reuse may lead to alterations in performance and cross-contamination risks. Dispose of the device after each single use; do not reuse. Avoid exposure to high temperatures and ultraviolet rays during storage. Adopt adequate precautions for elimination and disposal of the device and comply with the provisions of the laws in force on biologically hazardous waste.

STERILE - Sterilized by ethylene oxide

MATERIALS USED

Biocompatible and haemocompatible silicone.
LATEX-AND PVC-FREE MEDICAL DEVICE.

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Redax S.p.A.