





















## 18 872 B SysToe – Quick guide



















<i>Modif.</i>	<i>Rév.</i>	<i>Date</i>			
<b>04 090</b>	A	7/02/2011	First release		
<b>04 335</b>	B	28/08/2012	Release of new software version 23a		
<b>APPROBATION</b>					
<b>Quality</b>	<b>R&amp;D Atys</b>	<b>Sales</b>			

## Normal toe: MODE AUTO → two cuffs (occlusion & sensor cuff)


<b>1</b>	The supine patient should be resting for 10 min in a room between 22°C and 25 °C. The patient's feet should be warm (≈27°C)			
<b>2</b>	Wrap the occlusion cuff at the base of the toe. The cuff must not be too tight to prevent any residual pressure.			
<b>3</b>	It is compulsory to stick a piece of double-sided ring tape on the sensor. The tape must not cover the infrared cells (white rectangles).			
<b>4</b>	Stick the sensor on the pulp extremity of the toe and secure the sensor cuff. The sensor wire and the cuff tube go downwards as shown on the picture.			
<b>5</b>	Switch on the SysToe:  Check that the <b>AUTOMATIC MODE</b> is selected. If it is not the case, press on <b>MODE</b> .			
<b>6</b>	Press on <b>BRA.P</b> in order to input the arm pressure.			Input the pressure with the keyboard. Press on <b>RET..</b>
<b>7</b>	The displayed signal is pulsed or flat. Press on <b>START</b> . The measurement is performed automatically.			
<b>8</b>	Press on <b>STOP</b> when there is a clear and confirmed increase of the sensor signal. This action stops the examination before the complete deflation of the occlusion cuff.			Clear & confirmed increase of sensor signal
<b>9</b>	The opposite screen is displayed. Check the position of the vertical cursor. <b>It must be at the foot of the upslope of the sensor signal. If it is the case, press on VALID</b>			<b>If it is not the case</b> , move the cursor with the horizontal arrows. Press on <b>VALID</b> .
<b>10</b>	The systolic pressure and TBI are displayed.			<b>VALID</b> to save the measurement in the internal memory <b>BACK</b> Return to the previous screen. <b>NEW</b> Return to Step 7.

## Short toe: **MODE SEMI-AUTO** → only one cuff (occlusion cuff)

The SysToe is fitted with a short toe module. This module should be used only if the toe is too short to accommodate both the occlusion cuff and the sensor cuff.

<p>1</p>	<p>The supine patient. Wrap the occlusion cuff at the base of the toe. The cuff must not be too tight to prevent any residual pressure</p>			
<p>2</p>	<p>Remove the sensor from its cuff and stick a piece of double-sided ring tape on it.</p>			
<p>3</p>	<p>Stick the sensor on the toe pulp. The sensor wire and the cuff tube go downwards as shown on the picture.</p>			
<p>4</p>	<p>Switch on the SysToe:  Press on <b>MODE</b> in order to select the <b>SEMI-AUTOMATIC MODE</b>.</p>			
<p>5</p>	<p>Press on <b>BRA.P</b> in order to input the arm pressure.</p>			<p>Input the pressure with the keyboard. Press on <b>RET..</b></p>
<p>6</p>	<p>The displayed signal is pulsed or flat.</p>			
<p>7</p>	<p>Press on <b>START</b> The opposite screen is displayed</p>			
<p>8</p>	<p>Then press immediately strongly on the sensor as indicated on the opposite picture. Maintain the pressure as long as « <b>Keep press Until BIP</b> » is displayed.</p>			
<p>9</p>	<p>When this message disappears, release the pressure. The measurement is performed automatically. For the next steps, refer to the normal toe guide (step 8).</p>			

## CURVE ANALYSIS

	<p>The sensor detects the <b>return of blood flow</b> when the occlusion cuff is slowly deflated.</p> <p>When the blood flow comes back, the sensor signal curve shows a <b>significant change of slope</b> (the PPG curve clearly goes up). The systolic pressure is measured where the change of slope happens.</p> <p>At the end of the examination, a <b>vertical cursor</b> is automatically positioned on the blood flow return and the <b>systolic pressure (SYS)</b> is simultaneously displayed.</p> <p>NB: In case of positioning error, the vertical cursor can be adjusted by the user with the horizontal arrows.</p>
---	--

### Toe systolic pressure (SYS) and Toe Brachial Index (TBI)

#### Diagnosis of peripheral Arterial Disease

When arm pressure is input by the user, the SysToe calculates the TBI value.

$$\text{TBI} = \frac{\text{Toe systolic pressure (mmHg)}}{\text{Arm pressure (mmHg)}}$$

- TBI < 0.65 → PAD
- TBI > 0.65 → No PAD

#### Diagnosis of critical ischemia : SYS < 30 mmHg

#### Diagnosis of arterio venous hemodialysis access-induced hand ischemia

SYS < 60 mmHg or TBI < 0.4 are highly associated with hand ischemia.

**When the measured value is lower than the normal value, it is advised to perform a second measurement and even a third one to confirm the diagnosis.**



This short guide doesn't exempt you from reading the SysToe user's manual