REDAX[®]

IT

CATETERE TORACICO

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DESCRIZIONE

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INSTRUCTIONS FOR USE

DESCRIPTION

The Thoracic Catheter is a sterile, single-use, device. It is available in both the straight and the right-angled version, in silicone or PVC. All versions are available in different sizes. The PVC version has a specially shaped connector to facilitate inserting the catheter using aseptic technique. The proximal end of the silicone version is flute-lipped for easier insertion. Each catheter features a radiopaque line to allow radiographic visualisation. Depth markings, graduated from the last proximal eye, are provided at intervals to aid positioning and removal of the catheter. The product is supplied in sterile packaging.

INDICATIONS

The catheter is intended for use in medical situations that require drainage of air and/or fluid from the pleural and/or mediastinal space in the chest. Thoracic catheters are used for symptomatic patients whose pleural space has been invaded by air, fluid or pus; for patients with intrathoracic bleeding; and for post-thoracic surgery or heart surgery patients. Thoracic catheters are suitable for use as accessory devices for attachment to a thoracic drainage system for the collection and drainage of fluid and/or air from the pleural or mediastinal cavities following injury or surgical procedures.

Note: This device can only be used in safe conditions in the field of applications and using the procedures indicated on this operating instruction sheet for the specific type of product. The manufacturer declines all liability relating to improper use or use different to that indicated.

SUGGESTIONS FOR USE

Positioning the catheter.

- Using an aseptic technique, carefully remove the sterile thoracic catheter from its protective wrapping.
- An additional stab wound is made apart from the major incision.
- The catheter is inserted into the pleural space via the chest opening, and the proximal end extracted through the second stab wound in the chest wall.
- After properly positioning the catheter, make the connection to an underwater device, a suction source or a thoracic drainage device. In the PVC version, cut the blue cap on the end of the catheter (suitably disposing of the cut section). For flute-lipped silicone catheters, insertion is made through the proximal end that can be subsequently used for connection to the fluid collection system.
- It is suggested to use tape to secure the attachment of the connector to the drainage tube and also to the catheter.
- A suitable surgical technique should be used to secure the catheter to the chest wall of the patient.

Removing the catheter.

- Remove the tape or any other fastening and clamp the catheter to prevent air from penetrating the chest cavity.
- Deactivate the vacuum.
- Remove the catheter. Air leakage must be avoided during this process and a wound dressing must be immediately applied.

GENERAL WARNINGS

Do not use if the package has been opened or damaged. Carefully read this operating instruction sheet before using the product. Single-use product. Reuse may lead to alteration of performance and risks of cross-contamination.

CONTRAINDICATIONS

There are no known contraindications.

STERILE - Ethylene oxide gas sterilized

The product is sterile if the package is intact. Discharge after each single use, do not re-use. Avoid exposure to high temperatures and ultraviolet rays during storage. For discharging and disposal of the device, adopt suitable precautions and act in accordance with the provisions of the law in force for biologically hazardous waste

MATERIALS USED

PVC and Silicone biocompatible and hemocompatible.

LATEX-FREE MEDICAL DEVICE.

WARNINGS / PRECAUTIONS

- The right-angled catheter is designed for positioning between the chest wall and the lung; the 90° angle prevents the catheter from resting against the lung.**
- Before inserting the thoracic catheter, identify the introduction point in the chest wall and select the appropriate catheter size. Air drainage catheters are generally placed at the top front, while fluid drainage catheters are placed at the bottom rear.**
- Selection of the size of the device is based on physician preference and on the expected drainage quantity.**
- During removal, air leakage must be avoided and a wound dressing must be immediately applied.**
- Do not manipulate the device in any way nor make additional holes with the aid of cutting implements.**
- If fixing with a stitch, carefully avoid stitching the tube itself or cutting it with needles or sharp edges. Avoid a too tight stitch**

that might cause the drainage to rip.

- After positioning the drainage, ensure that the perforated section is completely inside the wound and that the connections are perfectly sealed.**
- The drainage must be removed by hand, pulling slightly and avoiding sudden movements. Avoid the use of metal instruments such as forceps or devices that might cause rupture of the drainage.**
- If the drainage is left in the site for very long periods of time, it might be difficult to remove it. Greater attention must therefore be paid during removal.**
- After removal ensure complete integrity of the drainage, an essential condition to exclude that drainage fragments caused by accidental rupture remain in the cavity.**
- It is recommended not to use the drainage tube for a period of more than 29 days.**
- Do not use roller clamps or similar metallic milking device on silicone tubes.**
- Do not use any abrasive tools, such as cotton flock or gauzes. Milking the silicone tubes by means of those tools can make damages on drain surface reducing the mechanical resistance.**
- Do not use any alcoholic or aggressive liquid to lubricate the surface of the tubes before milking. Redax strongly suggests Vaseline oil or cream to help the tube stripping or any other non-aggressive lubricant.**

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