

DRENTECH™ SIMPLY

Hybrid single bottle
drainage system



REDAX™
new drainage technologies born of experience

- Fast set-up aiming an easier drain management
- Reduce risk of cross-contamination
- Suitable for adult paediatric use

All-time clear visibility for air leakage




DRENTECH™ SIMPLY

looks alike the system


FAST LOCK-IN
CONNECTOR

VENTIL BALL-VALVE

INTEGRATED
UNDERWATER
SEAL VALVE



TRADITIONAL Single-bottle system



CHANGE YOUR HABIT

No more bottle replacement thanks to a capacious collection chamber (up to 1200 ml) with separated underwater seal valve.

DRENTECH™ SIMPLY & DRENTECH™ mobile REDLINE

- Simply and Redline are meant to be together, exalting each other features with an innovative design. Complete freedom for the Patient, independence from suction sources.
- Perfect choice for FAST TRACK surgery and early patient mobilization.



DrenTech Simply_1200 ml capacity _____code **10181**
Packaging: 12 units for box

DrenTech Simply_replacement _____code **10171**
Packaging: 12 units for box

DrenTech mini Simply_400 ml capacity _____code **10182**
Packaging: 12 units for box

DrenTech Mobile Redline
Kit portable chest drainage vacuum unit _____code **10188**
Packaging: 1 unit for box

Chest unit protection filter _____code **10539**
Packaging: 36 units for box



REDAX™

REDAX - Via Galileo Galilei, 18
46025 Poggio Rusco (MN) ITALY
Tel. +39 0386 830582 Fax +39 0386 51898
E-mail: info@redax.it - www.redax.it