

Suction Products, Suction Wands, Rigid Tip with No Block System with Plastic Handle, 24 Fr, S352-80



Livanova Rigid- and Flex-Tip Suction Wands are intended for intra-cardiac use. They are made of metal suction tubes with plastic handles. The Low Profile Suction Wand has a malleable low-profile metal suction tube with a plastic handle, pre-attached tubing and clamp, and terminates with a slip-on 1/4" connector. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

ORDER CODE	S352-80
SIZE	24 Fr
BARBED CONNECTOR	1/4"
AVAILABLE COATED	Depending on country registration
QTY PER BOX	10